



SkinPen®

Uncover the
true you
with SkinPen®

PREP. CARE. GLOW.

Preparing Your Skin for Renewal

Consultation & Skin Assessment

- ✓ Review your medical history, medications, and recent treatments, and inform your provider of any skin irritation in the treatment area.

Pre-Treatment Reminders

- ✓ Discontinue AHAs, Retinol, or Vitamin A.
- ✓ Avoid sun exposure and self-tanning products.
- ✓ If prone to cold sores, ask your provider about antiviral medication.

Provider Recommendations:

Supporting Your Skin's Recovery

Follow these steps for the best healing and results

Immediately After Treatment:

- ✓ Only use approved post-microneedling products.
- ✓ Apply Lift HG as needed to prevent drying out.
- ✓ Avoid makeup and sunscreen for 24 hours.
- ✓ No strenuous exercise or sweating for 72 hours.
- ✓ Stay out of direct sun for 72 hours.
- ✓ Use a SkinPen Biocellulose Masque to calm the skin.

In the Following Days:

- ✓ Skin may be pink and inflamed for 24-48 hours.
- ✓ Stay hydrated and avoid excessive heat (saunas, hot tubs).
- ✓ Apply RESCUE Calming Complex after 24 hours for soothing and skin support.

Maintaining Your Results

Use proper skincare to create an optimal environment for collagen remodeling, which can continue for up to six months post-treatment.

Provider Recommendations:

Maximizing Your Skin Rejuvenation

Clinical studies show that SkinPen® delivers optimal results with a recommended treatment series of 4 to 6 sessions, spaced 3 to 4 weeks apart.

Your scheduled SkinPen® treatment:

Thank you for choosing SkinPen®! Follow all post-care instructions for optimal healing.

Provider: _____

Phone: _____ Email: _____

SkinPen®

© Crown Aesthetics. All rights reserved.
SkinPen® is a registered trademark
of Crown Aesthetics. | F5SP485